

CORNING

Rock the Science of 3D! Corning Life Sciences Exhibitor Tutorial

Thomas Villani, Ph.D.

Chief Science Officer and Co-Founder
Visikol Inc.

Ann E. Rossi, Ph.D.


Applications Lab Manager
Corning Life Sciences

SLAS | February 2018






Corning Life Sciences:

An Innovator in Laboratory Products for Discovery Research




- Business Overview:
 - Market leader for over 100 years
 - Significant R&D investment for continuous innovation
 - Global manufacturing and distribution
 - Sales representatives worldwide including Technology and Field Applications Specialists
- Laboratory Products and Solutions for:
 - Cell Culture and Bioprocess
 - Drug Discovery
 - ADME/Tox
 - Genomics
 - Chemistry
 - Microbiology
 - General Laboratory Products



3D Image-Based Characterization of 3D Cell Culture Models Generated Using Spheroid Microplates


Tom Villani, PhD
Chief Science Officer
Visikol Inc.

Who is Visikol?




In Vitro Screening

- 3D Models
- High content confocal imaging
- Drug discovery assays



Whole Mount Analysis

- Vasculature mapping
- Whole brain imaging
- Whole mount IF, FP, ISH



Digital Pathology

- Slide scanning
- Quantitative analysis
- Feature extraction

- Biotechnology company that spun out of Rutgers University in 2016.
- Contract research organization focused on quantitative and 3D tissue imaging.
- Sells suite of patented Visikol[®] HISTO[™] reagents for use in research.



Tutorial Outline

Introduction to 3D Cell Culture

Corning® Ultra Low Attachment Spheroid Microplates


Background on Tissue Clearing

Limitations of High Content Screening


HepG2 Spheroid Application

Cells in 3D Culture More Closely Mimic *In Vivo* Physiology

Collagen-coated glass (2D)




Collagen gel (3D)



J Cell Sci 2012;125(Pt 13):3015-24.


LNCaP

C4-2



Adv Healthc Mater 2012;1(5): 590–99.

Cell Morphology
Structural Complexity
Phenotype



Pharmacol Ther 2016;163:94-108.

Spheroids, Microtissues, and Organoids

3D Cell Cultures

Spheroids

- Broad term for 3D agglomerations of cells, typically one or two cell types

Microtissues


- 3D cell cultures with multiple cell types designed to mimic the functionality of *in vivo* tissues

Organoids

- Highly complex organizations of cells designed to function similarly to *in vivo* organs

Common Techniques for Generating 3D Cultures


A. Matrix on-top and matrix- embedded



E. Ultra low attachment plates




B. Mat




Pharmacol Ther 2016;163:94-108.

Corning® Matrigel® Matrix
Corning® Ultra Low Attachment Surface



3D Cell Culture: Not a One-Size-Fits-All Solution



Realizing the Benefits of 3D Cell Culture can be Challenging


3D Strengths	3D Weaknesses
More physiological than 2D <ul style="list-style-type: none">• Cell morphology• Cell phenotype• More predictive of in vivo efficacy/toxicity	Culture conditions <ul style="list-style-type: none">• Variability, non-reproducible conditions with biological matrices• Hanging drop not amenable to high throughput without automation
Heterogenous cell populations	→ Culture and maintaining uniform size
Customizable environment	• Transfer necessary for assay, depending upon method
Growing body of work on 3D models, especially in oncology and toxicology	Assay methods/parameters not as well developed as 2D
Reagents specific for 3D	→ Imaging can be limited by 3D nature of samples
	Scalability

Corning® Spheroid Microplates


Corning Ultra-Low Attachment surface and unique round well-bottom design enable the formation and growth of a single, uniform spheroid per well with reproducible size.

Standard ANSI/SBS
footprint dimensions for
96- and 384-well formats

Black sidewalls to reduce
cross-talk and background
noise in fluorescent- and
luminescent-based assays



Clear bottom for
visualization and imaging




Corning Spheroid Microplate User Guide: Corning document CLS-AN-235


Reproducible and Routine Culture of Uniform Spheroids in Corning Spheroid Microplates

- Timing, seeding conditions dependent upon cell type, spheroid size

HT-29 Spheroids



Multicellular Spheroids




Scale bar = 1 mm

Traditional Methods for Interrogation


	Cost Per Well	Throughput	Limitations
Dissolution Assays	\$	High	Loss of spatial information
Chemical Analysis	\$	High	No spatial or cellular information
Widefield Imaging	\$\$	High	Only characterizes periphery of models
Histopathology	\$\$\$	Low	Tedious and low throughput
Confocal Microscopy	\$\$\$	Medium	Imaging depth limited to a few cell layers

High Content Screening (HCS)

- Phenotypic screening conducted on cells in well-plate format
- State of art has been to acquire single colorimetric readings, images or confocal z-stacks




U2OS actin, tubulin at 40X



Stitched image of neuron at 20X

Limitations in Imaging 3D Cell Cultures

- Opacity of 3D cell-cultures causes light scattering which limits imaging to $\sim 50\ \mu\text{m}$ deep, i.e. first 2-3 layers of cells
- Causes “eclipse” images common in microtissue confocal images
- Imaging bias: Majority of cells surveyed are the outer layer – most proliferative, most exposed to drug, oxygen, nutrients
- Histological sectioning of microtissues is cumbersome, not high throughput
- Many imaging endpoints currently assay total microtissue, unable to differentiate between cell populations
 - Diameter, total fluorescence, ATP-luminescence assays
- Myth: microtissue centers cannot be imaged because labels cannot penetrate to center



Tissue Clearing

- Numerous methods for clearing tissues have been developed over the last several years


Type of Technique	Technique	Refractive Index	Processing time for spheroids	Immuno-labeling Compatibility	Fluorescent Protein Compatibility	Tissue Integrity	Compatibility with well-plates	Clearing is reversible
Protein Hyper-hydration	Scale	1.38	Days	NO	YES	Expansion, very fragile	YES	NO
	ClearT2	1.44	Hours	YES	YES	Fragile	YES	NO
	CUBIC	1.38 - 1.48	Hours	YES	YES	Expansion, very fragile	YES	NO
Solvent Based	BABB	1.55	Minutes	YES	NO	Shrinkage	NO	NO
	i3DISCO	1.56	Minutes	YES	NO	Shrinkage	NO	NO
	Visikol® HISTO-M™	1.51	Minutes	YES	YES	No change	YES	YES
Aqueous RI Matching	FocusClear™	1.47	Hours	YES	YES	No change	NO	NO
	SeeDB/FRUIT	1.48	Hours	NO	YES	No change	YES	NO
Hydrogel Embedding	PACT/PARS	1.38 - 1.48	Hours/Days	YES	YES	Expansion	NO	NO
	CLARITY	1.45	Hours/Days	YES	YES	Expansion	NO	NO

Limitations of Known Clearing Techniques


- CLARITY + Hydrogel embedding techniques
 - Slow: requires days to weeks to clear even small tissues
 - Requires complicated processing technique – not amenable to medium-to-high-throughput
 - Not easily amenable to well-plate format
- Protein hyperhydration (CUBIC, Scale, SeeDB)
 - Slow: require days to weeks for clearing
 - Hyperhydration of proteins severely reduces tissue integrity
 - Not compatible with immunolabeling techniques
- Organic solvent-based (BABB, iDISCO, 3DISCO)
 - Use of highly lipophilic organic solvents not compatible with plastic well-plates

Limitations of Imaging Spheroids

Wide Field Microscopy



Non-cleared
57 cells detected




Cleared
151 cell detected


Confocal Microscopy

Non-cleared

Cleared



Clearing Allows Detection of Every Cell



Practical Considerations for Assaying Spheroids


- Generation of Spheroids
 - Choice of cell line
 - Plate choice
- Labeling techniques
- Considerations for High Content Imaging
 - Localizing spheroids
 - LED vs. Laser
 - Air objective vs. immersion objective
 - Confocal vs. widefield
 - Z-projections vs. Z-stacks

Labeling Techniques for Cleared Tissues

- Fluorescent probes
 - Work very well, very rapid in most cases
- Immunolabeling
 - Requires permeabilization to allow antibody penetration
 - Can label > 2 mm deep into brain tissue
 - Very powerful: label virtually any target in wild-type or genetically modified models
 - 1000s of antibodies available from vendors qualified for IF/IHC
- FISH
 - Requires protease pre-processing
 - Powerful tool to characterize gene expression
 - Expensive, often requires customized probes
- Fluorescent proteins
 - Excellent tool for interrogating gene expression, localization of target protein, etc.
 - Clearing provides excellent way to visualize fluorescent protein distribution

Considerations for Imaging: Localization

- Localization within the well
 - Flat-bottom plates cause spheroids to have random placement
 - Scan of entire well: ~72 hr. per 96-well plate
- Corning plates localize tissue to center of wells, drastically improving image throughput
 - Scan of entire plate: ~12 hr. per 96-well plate
 - Nearly 6x increase from U-bottom plates localizing tissue!



Considerations for Imaging: LED vs. Laser

- High content imaging systems utilize either LED or laser to excite fluorescent dyes
- Perkin Elmer Opera/Operetta, GE InCell 6000/6500, and Thermo CX7 LZR utilize laser illumination
- Molecular Devices, Thermo CX5/CX7 utilize LED
- Laser: better resolution, better penetration of light into tissue, no “foreshadowing” – slower, more expensive to maintain
- LED: lower resolution, causes “foreshadowing” – faster, less maintenance

Considerations for Imaging: Air Objective vs. Immersion

- For imaging spheroids, air objectives are perfectly compatible
- 10x air objective gives field of view (FOV) requiring ~36 tiles to cover well
 - Higher magnification requires even more scanning
- Immersion objective with high numerical aperture (NA) not required for routine spheroid work
- Myth: immersion objectives allow visualization of center of non-cleared spheroids
- Fact: spheroids must be cleared to visualize interior cells

Considerations for Imaging: Confocal vs. Widefield

- Confocal imaging limits out-of-focus light, resulting in sharper images with very narrow (~5-10 μm) depth of field
 - Volume images obtained by taking z-slice at various depth
 - With clearing, allows for interrogation of entire cell population
 - Much slower than widefield
 - Much crisper, sharper images than widefield
 - Smaller pinhole = sharper resolution = less light hitting detector = lower intensity image
- Widefield imaging
 - No depth information
 - With clearing, cross section of cell population can be observed
 - Much faster than confocal
 - Poorer image quality
 - Far more light hits detector, can detect weaker signals

Considerations for Imaging: Z-projections vs. Stacks

- Z-projections are a merge of images in a z-stack
 - Far smaller file size (2MB vs 50MB)
 - Loss of z-spatial information
 - Easier and much faster to analyze
 - Poorer quality than z-stack, better than widefield
- Z-stacks
 - Full set of coordinates for every detected cell
 - One image for each plane throughout spheroid
 - Very high image quality
 - Far more light hits detector: can detect weaker signals
 - Enables spatial response profiling

Data Processing: Open Source Software

- ImageJ: Open source image processing platform built on Java, from NIH
 - Clunky and difficult to learn, but **extremely powerful!**
 - Used for general processing (e.g. contrast adjustment, cropping, merging channels, etc)
 - Easy to process huge numbers of image files (>10k)
 - Easy to program using macros
- CellProfiler: Open source image processing platform for analyzing sets of images, from Broad Institute
 - Very powerful platform for cell-counting, colocalization, and other quantitative analyses of image sets obtained from high content imaging

Specific Application: HepG2 Cell Proliferation

- HepG2 spheroids were seeded at 1000 cells/well in Corning ULA U-bottom 96-well plates
- Grown for 2 days prior to dosing. Dosing day = day 0
- Dosed with paclitaxel at 1 μ M, 500 nM, 100 nM, 10 nM, 1 nM, as well as one row for vehicle control at day 0 and 2
- On day 4, spheroids were fixed with 10% neutral buffered formalin for 5 minutes

HepG2 Cell Proliferation: Processing

Pre-treatment (20 min.)

- Treatment with methanol/DMSO
- Treatment with PBS/Triton X-100



Labeling (1 hr.)

- Blocking
- Primary antibody (mouse anti-human Ki67)
- Secondary antibody (goat anti-mouse AlexaFluor488)



Clearing (5 min.)

- Dehydration
- Clearing with Visikol

HepG2 Cell Proliferation: Dose Response



HepG2 Cell Proliferation: Spatial Dose Response



Questions?

Poster #1036-B Today 5-6PM

Visikol Exhibit #1829

Corning Technical Program Handout

Visit Corning at Exhibit #1029 and
Rock the Science of 3D!



CORNING